



For reliable steaming / evaporation

equitron®

CONCENTRIC RINGS WATER BATH

Digital

Application:

- ❖ **EQUITRON** Concentric Rings Water Baths are ideal for all applications requiring gentle steaming work with Erlenmeyer flasks, glasses etc.

EQUITRON Advantages:

- ❖ Fitted with a water level indicating device.
- ❖ Optional constant water level device.
- ❖ A removable flat drop on lid for concentric rings enables the bath to be used for other work also.

Construction:

- ❖ Double walled, insulated, full stainless steel construction.
- ❖ Internal chamber, external chamber, drop on lid with holes of 75mm Ø all of stainless steel.
- ❖ Concentric rings in steps of 25 and 50 mm ID moulded from industrial grade non-corrosive material.

Standard Features:

- ❖ A precise Digital Temperature Indicator cum Controller having resolution and control accuracy of $\pm 1^\circ\text{C}$ in digital model.
- ❖ Temperature range: Ambient $+5^\circ\text{C}$ to 100°C .
- ❖ Industrial grade long life special heaters for efficient transfer of heat.



Model	Control	Capacity	Inner Dim. W x D x H (mm)	Concentric Rings	Heater	Packing Dimensions (WxDxH) / Gr. Weight
#6806DI	Digital	6 L	300 x 200 x 100	6 holes of 75 mm Ø with steps of 25 mm	800 W x 1	640 x 370 x 290 mm / 6 kg
#6812DI	Digital	12 L	400 x 300 x 100	12 holes of 75 mm Ø with steps of 25mm	800 W x 2	750 x 490 x 290 mm / 10.2 kg
#6824DI	Digital	24 L	600 x 400 x 100	24 holes of 75 mm Ø with steps of 25mm	1500 W x 2	990 x 580 x 280 mm / 15 kg

All electrical appliances work on 230 VAC, 50Hz, single phase.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series • Dry Bath • Colony Counter • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet • Test Tube Caps • Spreader • Cork Borer



MEDICA INSTRUMENT MFG CO

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West) Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ medica@medicainstrument.com • www.medicainstrument.com



0010_1005-3-1