



For convenience and economy

equitron®

ROUND WATER / OIL BATH IMMERSION SERIES Digital

Applications:

- ❖ *EQUITRON* Round Baths - with Digital control are ideal for reactions, extractions, analysis and tests.

EQUITRON Advantages:

- ❖ Double walled construction with efficient insulation system.
- ❖ Bath chamber of stainless steel.
- ❖ Outer body and control panel of MS duly finished in heat cured epoxy coating.
- ❖ On/off switch for mains and heater indication - on the front panel.
- ❖ Stainless Steel immersion heater - peripheral shape to eliminate direct contact with flask - even if flask touches the bottom.



Choice of Control:

- ❖ Digital Range: Ambient +10 to 250°C.
- ❖ Control accuracy: $\pm 2^\circ\text{C}$ or better.

Model	Round Flask Capacity	Bath size Ø x deep in mm	Temperature Range Water cum Oil - DWO	Heater	Packing Dimensions (WxDxH) / Gr. Weight
#8481ISDWO	1 lit.	Ø 165 x 90	Ambient + 10°C to 250°C	400 W	440 x 400 x 330 mm / 5.4 kg
#8483ISDWO	3 lit.	Ø 230 x 130		750 W	440 x 400 x 330 mm / 7 kg
#8485ISDWO	5 lit.	Ø 260 x 150		1000 W	550 x 460 x 370 mm / 8.1 kg
#8490ISDWO	10 lit.	Ø 350 x 200		1300 W	610 x 500 x 410 mm / 12 kg
#8492ISDWO	20 lit.	Ø 450 x 300		3000 W	770 x 640 x 480 mm / 19.4 kg

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series • Dry Bath • Colony Counter • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet • Test Tube Caps • Spreader • Cork Borer



Medica Instrument Mfg Co

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ sales@medicainstrument.com • www.medicainstrument.com



0011_1407-3-2