



equitron®

SEROLOGICAL WATER BATH

THERMOSTATIC

EQUITRON Advantages:

- ❖ Fully stainless steel, including bottom, makes it rust proof.
- ❖ Easy to set thermostat with 37°C calibrated mark.
- ❖ Ready to use - supplied with thermometer.

Construction:

- ❖ Double walled full stainless steel construction including bottom.
- ❖ Glasswool insulation for energy efficiency.

Standard Features:

- ❖ Capillary thermostat with $\pm 0.5^\circ\text{C}$ accuracy at 37°C.
- ❖ Energy efficient tubular heater.
- ❖ Range: Ambient plus 5°C to 65°C.
- ❖ Supplied with Thermometer, Clip and Mains Cord.



#6701TH

Model	Capacity	Inner Dimensions W x D x H (mm)	Racks (optional)		Optional SS Gabled Lid Model	Optional Float Ball Packs Model	Heater	Packing Dimensions (WxDxH) / Gr. Weight
			Kahn 30 #6205	Kahn 48 #6209				
#6701TH	4.2 L	270 x 130 x 120	-	1	#6901	#6990-50x1	500 W	440 x 400 x 330 mm / 4.4 kg
#6721TH	9.7 L	300 x 270 x 120	1	3	#6921	#6990-50x2	750 W	440 x 400 x 330 mm / 6.5 kg

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series
• Media Preparator • Table Top Autoclaves: Normal / Class B • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump •
BOD Incubator • Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series •
Dry Bath • Colony Counter • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren Apparatus • Slide Staining Stand /
Rack / Trough • Slide Trays • Slide Cabinet • Test Tube Caps • Spreader • Cork Borer



Medica Instrument Mfg Co

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ sales@medicainstrument.com • www.medicainstrument.com



0015_1506-2-2