



Save your environment & increase efficiency with

equitron®

CHILLED WATER CIRCULATOR - 4L

Applications:

- ❖ Ideally suitable for chilling and condensing operations. Ideal companion for Rotary Vacuum Evaporator, Viscosity, Spectrophotometry, Water Baths upto 28 lit and for general lab requirements of chilling in closed loop systems.

EQUITRON Advantages:

- ❖ Environment friendly CFC free hermetically sealed compressor.
- ❖ Sleek & compact design with easy maneuverability.
- ❖ Can be placed under workbench.

Construction:

- ❖ Stainless steel base & top. Panels of steel - powder coated.
- ❖ Four castors for easy movement.
- ❖ Supplied with 1 met x 2 Silicone tubings (8 x 11 mm) with insulation.
- ❖ Outer size W/D/H (mm): 560 x 410 x 640.

Standard Features:

- ❖ Suitable for closed loop systems.
- ❖ In-built Digital Temperature Controller cum Indicator.
- ❖ Temperature range upto -10°C (at ambient 25°C).
- ❖ Cooling capacity 750 W at 15°C.
- ❖ Reservoir capacity of approximately 4 lit.
- ❖ Circulation pump of app 15 LPM at zero head.



8502.ROT.000

Model	Description	Cooling Capacity	Packing Dimensions (WxDxH) / Gr. Weight
#8502.ROT.000	Chilled Water Circulator - 4L	750 W at 15°C	770 x 620 x 830 mm / 59 kg

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request. Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series • Dry Bath • Colony Counter • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet • Test Tube Caps • Spreader • Cork Borer



MEDICA INSTRUMENT MFG CO

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ medica@medicainstrument.com • www.medicainstrument.com



0076_1506-0-1