

# equitron<sup>®</sup> Essential Accessories for ROT EVA

## Vacuum Pump - PTFE coated diaphragm

### 5 mbar



- PTFE coated diaphragm pump, manufactured with German technology
- Maximum air displacement capacity: 38 ltr/min (2.28m<sup>3</sup>/Hr)
- Manufactured with German Technology
- Diaphragm: PTFE Coated
- Supplied with Solvent Trap
- Ultimate vacuum: 05 mbar



Model	Product	Ultimate Vacuum	Max Air Displacement	Packing Dimensions (WxDxH) / Gr Weight
#9072.005	Vacuum Pump	5 mbar	38 LPM	520 x 370 x 320 mm / 14 kg
#9072.002	Vacuum Pump	2 mbar	20 LPM	600 x 400 x 400 mm / 24 kg

# equitron<sup>®</sup> Essential Accessories for EVATOR

## Vacuum Pump - PTFE coated diaphragm

### 35 mbar



- PTFE coated diaphragm pump, manufactured with German technology
- Maximum air displacement capacity: 38 ltr/min (2.28m<sup>3</sup>/Hr)
- Manufactured with German Technology
- Diaphragm: PTFE Coated
- Supplied with Solvent Trap
- Ultimate vacuum: 35 mbar
- Supplied with solvent trap



Model	Product	Ultimate Vacuum	Max Air Displacement	Packing Dimensions (WxDxH) / Gr Weight
#9072.035	Vacuum Pump	35 mbar	38 LPM	520 x 370 x 320 mm / 14 kg

# equitron<sup>®</sup> Vacuum Pump - Variable

## 2 mbar / 20 LPM



- PTFE coated variable drive energy efficient diaphragm pump manufactured with German technology
- Digital controller with display of actual and set values.
- Supplied with Solvent Trap



Model	Product	Ultimate Vacuum	Max Air Displacement	Packing Dimensions (WxDxH) / Gr Weight
#9072.002.VD	Vacuum Pump - Variable	2 mbar	20 LPM	910 x 640 x 510 mm / 30 kg

