



No Heater Burnouts

# equitron®

## WATER STILL

With Thermostatic cut-off

### Applications:

- ❖ **EQUITRON** Water Still is ideal for all laboratories where distill water free from dissolved salts, clear, odourless and tasteless water is required.

### EQUITRON Advantages:

- ❖ Thermostatic cut-off (self resetting type) for assured low water detection - with visual indication.
- ❖ Unique water level maintaining device improves efficiency.
- ❖ Compact & user friendly space saving design - occupies less space & can be table or wall mounted.
- ❖ Easily openable Boiling Chamber as well as Condenser for easy cleaning.

### Construction:

- ❖ Fully stainless steel.
- ❖ Press fit boiling chamber with internal gasket & condenser.

### Standard Accessories:

- ❖ Silicone tubing for distill water output + for flow of water from Condenser to Boiling Chamber.
- ❖ PVC tubing for waste water + raw water with tap adapter and a hose clip.
- ❖ Flow rate meter plus Hoffman clip for raw water flow monitoring.
- ❖ Mounting screws for Water Still and Flow meter.



Model	Capacity	Heater	Packing Dimensions (WxDxH) / Gr. Weight
#7655M	4 lit / hr	3.0kW	550 x 460 x 370 mm / 6.8 kg

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request. Due to continuous improvement in design, the product supplied may have modified features.

**Our Other Products:** Autoclave Vertical: Standard / PAD / Fully Automatic, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series • Dry Bath • Colony Counter • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet • Test Tube Caps • Spreader • Cork Borer.



**Medica Instrument Mfg Co**

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)  
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996  
✉ sales@medicainstrument.com • www.medicainstrument.com



0015\_1605-3-1