



EQUITRON[®]

Petri Dish Sterilizing Box

- ❖ Fully stainless steel construction.
- ❖ Petri dish box is seamless, die pressed.
- ❖ With removable stainless steel rack (dish carrier).



| Model | Product | Size | Accommodates Petri Dish |
|---------|-----------------------------|--------------|---------------------------------|
| #6620-H | Petri Box with carrier - SS | Ø7.5 x 19 cm | Ø60 ~ 65 mm Approx. 12 nos. |
| #6621-H | Petri Box with carrier - SS | Ø11 x 25 cm | Ø90 ~ 100 mm Approx. 12 nos. |

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • BSC Clean Bench • Water Still / Stone-fin / Automatic Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Shaking / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series / ECO Series • Ultrasonic Bath: Digital / LCD • Dry Bath • Colony Counter • Hot Plate with Magnetic Stirrer • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren & Wintrobe Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet



MEDICA INSTRUMENT MFG. CO.

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ sales@medicainstrument.com • www.medicainstrument.com



0147_1201-0-1