


equitron®

Slide Cabinet - Individual



- Expandable & Flexible storage system for slides
- Each stack able cabinet stores 2520 slides (individual) in just 1.25² ft
- Can stack upto 3 cabinets - one above the other
- Drawer does not fall down even when fully pulled out

Model	Product	capacity
#6005	Slide Cabinet - Individual Storage - 36 drawers	2520 slides
#6050	Base for above	Stack upto 3 cabinets

equitron®

Slide Trays / Stand



- Viewing tray: aluminium embossed and powder coated
- Drying tray: stainless steel glossy inish - drains off water easily
- Staining Stand (ghodi): full stainless steel rigid construction

Model	Product	capacity / description
#6105	Viewing tray - AL	24 slides
#6115	Drying tray - SS	4 grooves of 300 mm
#6121	Staining Stand (Ghodi) - SS	280 mm long

equitron®

Slide Staining Racks & Troughs



- Stainless steel die punched rack
- Stainless steel joint less staining trough with rounded edges for easy cleaning

Model	Product	capacity
#6142	Slide Staining Rack - SS	10 slides
#6145	Slide Staining Rack - SS	25 slides
#6152	Slide Staining Trough - SS	10 slides
#6155	Slide Staining Trough - SS	25 slides

equitron®

Petri & Pipette Sterilising Boxes (Cannisters)



- Petri Boxes in stainless steel joint less die pressed body & lid
- Petri Box supplied with internal Petri Dish Carrier suitable for approx. 10 dishes
- Pipette Box in full stainless steel construction
- Unique autoclavable tip saver cushion at bottom in Pipette Boxes



Model	Product	Dimensions
#6601	Pipette Box SS	Ø 60 x 380 mm (15 x 2½"Ø)
#6605	Pipette Box SS	Ø 50 x 450 mm (18 x 2" Ø)
#6620	Petri Dish Box for upto Ø 65 mm (outer) plates	Ø 75 x 190 mm (7½ x 3" Ø)
#6621	Petri Dish Box for upto Ø 100 mm (outer) plates	Ø 110 x 250 mm (10 x 4¼" Ø)

equitron®

Test Tube Racks Double Tier - SS



- Stainless steel die punched burr free holes
- Double Tier racks - Top & centre punched, bottom plain



Model	Product	Dimensions (hole conig)
#6205	Kahn - 30	Ø 13 x 30 holes (3x10)
#6209	Kahn - 48	Ø 13 x 48 holes (4x12)
#6235-12	2T Ø20x12	Ø 20 x 12 holes (2x6)
#6235-24	2T Ø20x24	Ø 20 x 24 holes (3x8))
#6235-48	2T Ø20x48	Ø 20 x 48 holes (4x12)
#6242-12	2T Ø26x12	Ø 26 x 12 holes (2x6)
#6242-24	2T Ø26x24	Ø 26 x 24 holes (3x8)
#6245-12	2T Ø32x12	Ø 32 x 12 holes (2x6)
#6248-12	2T Ø40x12	Ø 40 x 12 holes (2x6)
#6248-24	2T Ø40x24	Ø 40 x 24 holes (3x8)

equitron®

Test Tube Caps



- Stainless Steel Deep Drawn Test Tube Caps
- Test tube Cap for 38 mm tube extra long and heavy weight



Model	Product
#6390-27	Test Tube Cap SS - for 25 mm tube
#6390-3954	Test Tube Cap SS - for 38 mm tube



Medica Instrument Mfg Co
www.medicainstrument.com